

Glow discharge system for TEM grids

- Dual or single chamber
- Hydrophilic/hydrophobic and negative/positive modes
- Fully automatic, short process times
- Intuitive touch screen control
- Safe vapour delivery using septum-sealed vials
- Automatic valving between chambers to prevent cross-contamination
- Quick and easy sample loading
- Controlled venting to prevent sample disturbance
- Consistent, reliable results
- Three-year warranty



Rapid, reliable results with the GloQube™

The GloQube™ is a compact, easy to use, stand-alone glow discharge system.

The primary application of the GloQube™ is the hydrophilisation (wetting) of carbon-coated TEM support films and grids which otherwise have the tendency to be hydrophobic. Glow discharge treatment with air will make film surfaces negatively charged and hydrophilic and allow the easy spread of aqueous solutions. This and other processes are outlined below.


Glow discharge processes			
Surface state	Charge	Atmosphere	Typical applications
Hydrophilic	Negative	Air	Carbon coated TEM grids
Hydrophilic	Positive	Air – with magnesium acetate post-treatment	Nucleic acid adhesion to carbon films
Hydrophobic	Positive	Alkylamine	Proteins, antibodies and nucleic acids
Hydrophobic	Negative	Methanol	Positively charged protein molecules (e.g. ferritin, cytochrome c)

Single or dual chamber, safe handling of reagents

The GloQube™ is available in two formats:

The GloQube-S has a single “clean” chamber designed for applications requiring hydrophobic-hydrophilic conversion, typically using air as the process gas.

The GloQube-D has two independent vacuum chambers: a clean chamber and a vapour chamber, which is designed for use with reagents such as methanol and alkylamine. With operator safety firmly in mind, reusable septum-sealed reagent vials are used. Loading and removing reagents is convenient and reliable – the vial, located in its holder, is inserted into a shielded needle using a simple bayonet fitting.



Start-up screen



Clean chamber



Vapour chamber



Vapour delivery system




Easy sample loading, fast turnaround times

Each chamber can accommodate two 25 x 75 mm glass microscope slides. Loading could not be easier using draw-style chamber doors and specimen stages. The stages are height adjustable and fitted with removable glass slide holders. For additional convenience – and to allow easy access for chamber cleaning – the stages can be completely removed.




Touch screen control – rapid data input, simple operation


The intuitive touch screen allows multiple users to rapidly input and store preferred process “recipes”. Typical default glow discharge protocols are loaded as standard. Additionally, help files and useful maintenance data such as system on time and time since last clean are readily available to the operator. An Ethernet communications port is included for software updates.



Stored profiles



Selecting a new profile



A typical process run




GloQube™ clean chamber during a process run

Vacuum, automatic valving and controlled venting

The GloQube-D has automatic valving between chambers which maintains cleanliness by preventing cross-contamination. At the end of a process run, automatic soft venting to atmosphere through filtered inlets ensures TEM grids are not disturbed.

Both GloQube-S and GloQube-D require a single vacuum pump working in the 0.1 to 1 mbar range. A typical pump time to operational vacuum is 60 seconds.



Pumping cycle



GloQube-D and optional Pfeiffer DUO 6 rotary pump

Ordering information

- **025235** GloQube-D. Dual chamber glow discharge system. 025678 accessory kit, including: mains power lead, rotary pump power lead, 07803 oil mist filter and KF16 clamp, 750 mm long flexible stainless steel vacuum tube with clamps, fuses
10 x 025266 glass vials, 3 x vial caps and sealing washers, 1 x 025345 needle (spare)
Vacuum pump to be ordered separately
- **025261** GloQube-S. Single chamber glow discharge system. 025938 accessory kit, including: mains power lead, rotary pump power lead, 07803 oil mist filter and KF16 clamp, 750 mm long flexible stainless steel vacuum tube with clamps, fuses. Vacuum pump to be ordered separately

Vacuum pumping

- **13034** 5 m³/hr Pfeiffer DUO 6 two-stage rotary vacuum pump with 07803 oil mist filter
- **07803** Oil mist filter (spare)

Options, accessories and spares

- **026032** Single chamber upgrade kit. For retrospective conversion of GloQube-S to GloQube-D
- **025195** Microscope slide tray
- **025266** Glass vial (packet of 10)
- **025267** Glass vial caps (packet of three)
- **025345** Needle
- **024755** Door seal

Specifications

		GloQube-S Single chamber	GloQube-D Dual chamber
Power and processes			
Plasma current	1-40 mA	✓	✓
HV power supply	30 W	✓	✓
Maximum voltage	800 V	✓	✓
Electrode polarity – clean chamber	DC glow positive DC glow negative	✓	✓
Electrode polarity – vapour chamber	DC glow positive DC glow negative	NA	✓
Sample stage	125 x 100 mm (4.9" x 3.94") with location for two 25 x 75 mm (1" x 3") glass slides	x 1	x 2
Sample stage operational heights	Adjustable 12.5 mm (0.5"), 22.5 mm (0.9") or 35 mm (1.38")	✓	✓
Pump hold time	0-24 hours	✓	✓
Process time	1-600 seconds	✓	✓
Safety			
Chamber vent inlets	Filtered air inlets with slow vent to minimise sample disturbance	x 1	x 2
On-board reagent storage	Reagents (e.g. methanol or alkylamine) are contained in reusable sealed glass vials to minimise exposure to hazards. (GloQube-D only)	NA	✓
High voltage safety interlocks	Hardware safety interlocked and software for process control	✓	✓
Vacuum			
Vacuum control	Integrated pirani gauge	x 1	x 2
Working vacuum range	0.1 to 1 mbar	✓	✓
Vacuum pump minimum requirements	6 m ³ /hr, 3600 l/m, 0.03 mbar ultimate vacuum. Inlet flange: KF 16	✓	✓
Pumping time	Typical pump time to an operational vacuum of 0.27 mbar in 60 seconds	✓	✓
Vacuum isolation	Isolation valves to switch vacuum and prevent process chamber cross-contamination	x 1	x 2
User interface			
User interface	Full graphical interface with touch screen buttons and controls. In addition to displaying profiles, parameters, help screen and maintenance information are available	✓	✓
Profiles and profile logging	Capability to store 100 user profiles (name, date, time, vacuum, current and polarity)	✓	✓
Dimensions and communications			
Chamber size	100 mm W x 100 mm H x 127 mm D (3.94" x 3.94" x 5")	x 1	x 2
Instrument size	336 mm H x 364 mm D (13.2" x 14.3")	350 W (13.8")	366 W (14.4")
Instrument weight	19.5 kg (42.9 lbs.) (GloQube-D)	TBA	19.5
Pump (optional)	391 mm W x 127 mm D x 177 mm H (15.4" x 5" x 7")	✓	✓
Pump weight	16 kg (35.3 lbs)	✓	✓
Footprint with optional pump	366 mm W x 600 mm D x 336 mm H (14.4" x 23.6" x 13.2")	✓	✓
Power requirements	120 V 60 Hz, 15 A or 230 V 50 Hz, 10 A	✓	✓
Instrument power rating	100-240 V AC 60/50 Hz 700 VA including pump, IEC inlet	✓	✓
Optional pump power rating	115/230 V 60/50 Hz 450 W	✓	✓
Communication port	Ethernet port for instrument software updates	✓	✓



Quorum Technologies Ltd

Judges House, Lewes Road, Laughton, Lewes, East Sussex BN8 6BN UK

T: +44(0)1323 810981

E: sales@quorumtech.com www.quorumtech.com