

### **Uniblitz® FS25**

25mm Uni-Stable Optical Shutter

### Overview

The Uniblitz FS series shutters are designed and optimized to operate directly from +5VDC and do not require a separate driver. Removing the +5VDC (0VDC) closes these shutters. A simple control circuit can be used to operate these shutter devices from a TTL trigger pulse. This control can be also accomplished with our new VLM1 TTL control interface (available soon).

This low-cost innovation provides the reliability of Uniblitz shutters (typical lifetime >300K cycles) at a single operating voltage.

### **Need Support?** Please <u>visit our website</u> or email us at <u>info@uniblitz.com</u>. Tel: <u>585-385-5930</u> | Toll-Free: <u>800-828-6972</u> | Fax: <u>585-385-6004</u> | 803 Linden Ave. Rochester, NY 14625 Updated 1/20 | Datasheet Version 5.2 | ©2019 Vincent Associates

### **Key Features**

- 25mm aperture
- Default closed operation, +5VDC opens the shutter, 0VDC closes the shutter
- RoHS Compliant
- Transfer time on opening:
   9.0 milliseconds
- Total opening time:
  - 16.0 milliseconds

## **Product Options**

### FS25 2 3 4 5 6 - 7

#### Ex: FS25S2C0L-EC

impregnated plastic (C-PET)

• T: Low Energy (Teflon<sup>®</sup>)

• ZM: High Energy (AlMgF2)<sup>2</sup>

**2** Voltage:

4 Blades: 1

**6** Connector:

• L: 18" flying leads

• **C:** Black carbon


• **S:** Standard 5VDC

- 1 Shutter Series:
- FS25: Normally Closed
- FSR25: Normally Open
- **3** Housing:
- 1: Un-Housed
- 2: Half-Housed
- 3: Fully-Housed
- **5** Electronic Sync:
- 0: Omitted
- 1: Included
- **7** Encapsulated Coil:
- EC: Included <sup>3</sup>
- Leave blank if not required

<sup>1</sup> Various mounting methods are available depending on housing option – see Key Features on website

<sup>2</sup> Other blade coating options may be available by special order.
<sup>3</sup> Input side only; Teflon<sup>®</sup> coating is on opposite side to protect shutter blade surface. Light source must be input to the reflective side only.
<sup>4</sup> Vacuum compatibility up to 10E-10 Torr

# Shutter Timing



### FS25 (w/ 5VDC and C-PET blades)

#### Time (msec.)


| 0 - A | Delay time on opening after current applied | 7.0   |
|-------|---|-------|
| A - C | Transfer time on opening                    | 9.0   |
| 0 - C | Total opening time                          | 16.0  |
| C - E | Min. dwell time with min. input pulse       | 15.0  |
| B - F | Min. equivalent exp. time                   | 30.5  |
| E - G | Transfer time on closing                    | 22.0  |
| A - G | Total window time                           | 53.0  |
| MET   | Min. exposure time                          | 30.0  |
| TEP   | Typical exposure pulse                      | >30.0 |

## **Technical Specifications**

| Coil<br>Resistance | Voltage to Open | Hold Voltage |
|--------------------|-----------------|--------------|
| 12 Ω               | +5 VDC at 425mA | +5 VDC       |


<sup>1</sup> (Continuous/Burst) Continuous frequency rating specified at shutter's minimum exposure pulse. Burst frequency rating specified for four (4) seconds maximum with one (1) minute minimum between bursts.

| Series | <b>Weight</b>            | Operating    | Max. Opening | Max. Closing | Max. Freq. of          | Number of      |
|--------|--------------------------|--------------|--------------|--------------|------------------------|----------------|
|        | (Unhoused/Half/Housed)   | Temp.        | Bounce       | Bounce       | Operation <sup>1</sup> | Shutter Blades |
| FS25   | 22.0 g / 57.0 g / 92.5 g | -40 - +65 °C | 15%          | 15%          | 5 Hz / 10 Hz           | 5              |



The synchronization system for FS shutter devices incorporates a small magnet mounted to the driving mechanism and a Hall effect sensor. When the device achieves approximately 80% of full open, the magnet causes the Hall effect sensor to change state, producing a signal to indicate that the shutter has switched to the active state. Shown to the left is the FS series shutter schematic which incorporates the electronic synchronization system. **There is no connection to the designated synchronization pins when an electronic sync. is not selected.** 

# **Typical Control Circuit**




This simple control circuit can be used to operate the shutter device from a TTL trigger pulse.<sup>1</sup> This control can be also accomplished with our new VLM1 TTL control interface, which will be available soon.

<sup>1</sup> The shutter will not operate directly from a TTL signal.


### Technical Drawings - FS25 (Half-Housed - FS25S2C0L)

#### NOTES:

- 1. UNITS: [MM] INCH.
- 2. LIMITED HIDDEN LINES SHOWN FOR CLARITY.
- 3. LEADS WIRES REMOVED FOR CLARITY.
- 4 Shutter actual aperture is 25.4mm, in mount aperture is reduced by female C-mount to 0.965 [24.5mm].



### Technical Drawings - FS25 (Un-Housed - FS25S1C0L)



## Technical Drawings - FS25 (Housed - FS25S3C0L)

#### NOTES:

- 1. UNITS: [MM] INCH
- 2. LIMITED HIDDEN LINES SHOWN FOR CLARITY.
- 3. SHUTTER APERTURE IS 25.4MM REDUCED TO (Ø.965 [24.5MM]) WHEN INSTALLED INTO HOUSING.
- 4. OPTIONAL REFLECTIVE SURFACE OPPOSITE ACTUATOR COIL SIDE.

